

Ur24T TrueClr™

Introducing TrueClr Urinary Collection Apparatus (UCA), the world's first completely external catheter system designed to comfortably empty the bladder while eliminating the risk of a Catheter Associated Urinary Tract Infection (CAUTIs).

Innovatively Designed

- Comfortable and conforms to the human anatomy.
- Empties the bladder without pain or discomfort.
- Designed for both males and females.
- Preserve dignity, save time & manage costs.
- Easy to use, portable & reusable.

How Does TrueClr™ Work?

The Ur24T TrueClr external catheter operates by actively removing urine from the bladder through light suction. The TrueClr external urinary catheter suction is powered by an aspirator. Various models of aspirators can be used in an assisted living setting, or the catheter can be powered by a wall aspirator that is standard in most hospitals.

