

# DINAMIC OR<sup>+</sup>

Distributed by Maclean Health



## THE ULTIMATE PROCEDURAL SOLUTION

Finally, a Customizable Back Table Solution that will Optimize Your OR.

# COMPLEX PROBLEMS IN THE OR...



---

Inconsistent Setups Cause  
Errors and Waste Time



---

Staff Turnover



---

Onboarding & Knowledge  
Transfer



---

Setup & Break Down Time  
with Lack of Workspace



---

Tray & Instrument  
Inventory Issues



---


Increased Demand  
for OR Time

# ...REQUIRE A SIMPLE, COMPLETE SOLUTION.

+ EFFICIENCY

33% TIME  
REDUCTION


ADJUSTABLE  
BACK TABLE



+ SAFETY

88%  
DECREASED  
TRAY MOVEMENTS


DURABLE DRAPE




+ EDUCATION

PROVIDE JUST-IN-TIME  
**TRAINING**  
FOR SURGICAL TECHS

WORKFLOW  
MANAGEMENT APP




WATCH NOW: "THE DINAMICOR SUITE OF  
SOLUTIONS DRIVEN BY FORMULA 1"



# ADJUSTABLE BACK TABLE SOLUTION

With an ergonomic design and adjustable tiers, the DinamicOR Back Table caters to the preferences of your surgical team, minimizing strain and enhancing efficiency. With a similar footprint to your standard table, the vertical, multi-shelf design allows you to maximize the workspace in your OR. Configurable to different procedures, the 4-in-1 design allows for easy adjustment of shelf positions. Choose between two sizes, accommodating up to 8 or 12 trays.



## TRIPLE WIDE

### Full Table

Dimensions: 68.0\"L x 28.3\"D  
Height Range: 49\" - 55\"

### Main Surface

Dimensions: 68.0\"L x 26.3\"D  
Height from Floor: 30\" - 36\"  
Weight Capacity: 300lbs

### Angled Surfaces

Top Shelf: 68.0\"L x 11.3\"D  
Middle Shelf: 68.0\"L x 12.3\"D  
Weight Capacity: 150lbs

## DOUBLE WIDE

### Full Table

Dimensions: 45.5\"L x 28.3\"D  
Height Range: 49\" - 55\"

### Main Surface

Dimensions: 45.5\"L x 26.3\"D  
Height from Floor: 30\" - 36\"  
Weight Capacity: 300lbs


### Angled Surfaces

Top Shelf: 45.5\"L x 11.3\"D  
Middle Shelf: 45.5\"L x 12.3\"D  
Weight Capacity: 150lbs

## 4-TABLES-IN-1

Foldable Shelving Configurations






## UNIQUE FITTED DURABLE DRAPE

Thick, durable, and tear-resistant, the DinamicOR custom-fitted drape is uniquely designed for the back table shelving, that never requires 'double-draping'. Complete with numbered tray placement indicators for fast and efficient setup.



## REPEATABLE TRAY PLACEMENT

The durable Drape incorporates numbered placement icons that line up with the shelves for guided setup. This allows you to repeat your preferred placements every time.



## EASY TO INSTALL


The DinamicOR durable Drape is specifically designed to fit the adjustable Back Table and is tightly secured with customized wire tape. The durable Drape is the thickest drape on the market, which means you don't need to use multiple drapes per procedure.



# THE WORKFLOW MANAGEMENT APP

The DinamicOR+ Workflow Management App allows you to map your preferred tray layouts for any procedure and then repeat them as needed. This allows the surgical team to find the instruments they need exactly where and when they need them.

## TABLE SETUP MODE: PLACEMENT & VERIFICATION



### RISK MITIGATION

Verify all necessary equipment is in the room prior to starting the case.


### IMPROVE EFFICIENCY

The DinamicOR+ App's guided walkthroughs reduce setup time by up to 33 percent.

### REPEATABILITY

Set up the back table the same way every time no matter who is in the room.

## SURGICAL WORKFLOW MODE: ONBOARDING & EDUCATION



### SURGICAL CHECKLISTS

Provide staff just in time education during the case

### TIME SAVINGS

Guides staff through every step of setup. Complex and unfamiliar procedures can be mapped and planned to reduce OR errors.

### PROCEDURAL SEQUENCING

Plan every procedure to have the exact needed instruments every time, ensuring fast and reliable setup.

## EMPOWERING EDUCATION & EFFICIENCY

- Customize procedures for smoother workflows
- Share knowledge with documented workflows
- Simplify staff onboarding
- Ensure accurate instrument placement with a reference database
- Reduce risk by verifying instrument availability

# UNLOCKING OPERATIONAL EXCELLENCE

---

## EFFICIENCY

- » Reduce setup and breakdown time by up to 33%.
- » Maximize OR workspace with vertical multi-shelf design.
- » Enhance workflow with 92% instrument visibility.

---

## SAFETY

- » Ensure a sterile environment with a durable, single-drape solution.
- » Mitigate risks with the app's verification of all necessary instruments.
- » Reduce staff injuries with 6" adjustable table height.

---

## ONBOARDING & EDUCATION

- » Streamline onboarding with guided tray setups.
- » Facilitate knowledge transfer with custom procedural mapping.
- » Optimize education and training with a user-friendly Workflow Management App.

## OUR MISSION IS TO IMPROVE OR ORGANIZATION, SAFETY, AND EFFICIENCY

Our mission is to optimize operating room organization with an ergonomic and intuitive Workflow Management System – improving the efficiency and reproducibility of surgical processes which will enable perioperative staff to focus on providing thoughtful patient care.



### Watch the Table Setup Demo

More Info contact: Sheri Hottle

Call: +1 (317) 500-0392

Email: [shottle@macleanhealth.com](mailto:shottle@macleanhealth.com)

URL: [www.DinamicOR.com](http://www.DinamicOR.com)

SCHEDULE A LIVE DEMO TO LEARN MORE  
ABOUT OUR CUSTOMIZABLE SOLUTIONS

Distributed by Maclean Health

