



BATHROOM SAFETY



Three in One Commodes

The Lumex® three in one commodes are designed to be used as a bedside commode, raised toilet seat, or toilet safety frame. Each one comes with commode pail, lid, and splashguard. The backbar detaches so that commode frame can be used as a toilet safety frame.

Steel Three in One Commodes

- Strong, stable, durable welded-steel frame
- Powder coated frame is attractive and easy to clean
- Item No. 7103, 7103A-4



Folding Steel Three in One Commodes

- Folds for convenient storage and transport
- When folded, the arms can be locked into place with the back bar for easy transport
- Item No. 7108A, 7108R-1



Aluminum Three in One Commodes

- Lightweight, easy-to-transport aluminum frame
- Contoured armrests enhance user comfort
- Item No. 2195A-4, 2215B-2



Drop Arm Steel Three in One Commodes

- Drop arm is easy to release and operate by patients lacking manual dexterity
- Drop arms are independent and lock automatically in raised position
- Padded armrests provide comfort
- Item No. 6433A, 6433A-2



Item No.	Width Between Arms	Overall Width	Overall Depth	Seat Width	Seat Depth	Seat Height	Weight Capacity*	Retail Pckg
7103	18.5 in / 46.9 cm	24.5 in / 62.2 cm	19 in / 48.2 cm	14 in / 35.5 cm	16.5 in / 41.9 cm	17 in - 27 in 43.1 cm - 68.5 cm	300 lb / 136 kg	
7103A-4	17.5 in / 44.4 cm	24 in / 60.9 cm	19 in / 48.2 cm	14 in / 35.5 cm	15 in / 38 cm	15.5 in - 23.5 in 39.3 cm - 59.6 cm	300 lb / 136 kg	
7108A	20 in / 50.8 cm	24 in / 60.9 cm	22 in / 55.8 cm	14 in / 35.5 cm	16.5 in / 41.9 cm	17 in - 23 in 43.1 cm - 58.4 cm	300 lb / 136 kg	
7108R-1	18 in / 45.7 cm	22 in / 55.8 cm	22 in / 55.8 cm	14 in / 35.5 cm	16.5 in / 41.9 cm	17 in - 23 in 43.1 cm - 58.4 cm	300 lb / 136 kg	✓
2195A-4	18.5 in / 46.9 cm	26.5 in / 67.3 cm	18 in / 45.7 cm	13.5 in / 35.5 cm	14.5 in / 37 cm	17.5 in - 21.5 in 44.4 cm - 54.6 cm	300 lb / 136 kg	
2215B-2	20 in / 50.8 cm	24 in / 60.9 cm	21.5 in / 54.6 cm	13.8 in / 34.9 cm	16.5 in / 41.9 cm	21.5 in / 54.6 cm	300 lb / 136 kg	
6433A, 6433A-2	18 in / 45.7 cm	23.7 in / 60.3 cm	19 in / 48.2 cm	14 in / 35.5 cm	16.7 in / 43 cm	18.7 in - 22.8 in 47.6 cm - 57.7 cm	300 lb / 136 kg	

Specifications are subject to change without notice

*EVENLY DISTRIBUTED

Platinum Bath Seats

The Lumex® Bath seats are designed to provide comfort and stability to those needing to sit down in a shower or bathtub. These seats are lightweight, durable, and have a rust-resistant anodized aluminum frame. The seat height is adjustable in 1"

increments in order to best adapt to a user's height. The Lumex bath seats have one of the widest seats on the market for comfort and are available with or without a backrest.

Assembled Bath Seats

- Comes assembled with tool free back attachment
- Item No. 7921A-4, 7921A-1 (With Backrest)
- Item No. 7931A, 7931A-1 (Without Backrest)

Knock-Down Bath Seats

- Tool-free design with easy snap together parts
- Attractive multiple colors available
- Item No. 7921R, 7921RB, 7921RB-1, 7921KD-4, 7921KD-1 (With Back)
- Item No. 7931R, 7931RB, 7931KD-4, 7931KD-1 (Without Back)

Bath Seat With Arms

- Removable support arms provide for easy positioning
- Item No. 7921RA (With Backrest)
- Item No. 7931RA (Without Backrest)

Bath Seat Basket

- Keeps bath items accessible
- Helps avoid unnecessary standing and reaching
- Rust resistant
- Item No. 7921R-BSKT-1



Item No.	Width Between Arms	Overall Width	Overall Depth	Overall Height	Seat Width	Seat Depth	Seat Height Range	Back Support Height (Above Seat)	Weight Capacity*	Retail Pckg
7921A-4	—	20 in / 50.8 cm	18 in / 45.7 cm	28 in - 33 in 71.1 cm - 83.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	15 in - 20 in 38.1 cm - 50.8 cm	13 in / 33 cm	350 lb / 158.8 kg	
7931A	—	20 in / 50.8 cm	18 in / 45.7 cm	28 in - 33 in 71.1 cm - 83.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	15 in - 20 in 38.1 cm - 50.8 cm	—	350 lb / 158.8 kg	
7921KD-4, 7921KD-1	—	20 in / 50.8 cm	19 in / 48.2 cm	27 in - 33 in 68.5 cm - 83.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	14 in - 20 in 35.5 cm - 50.8 cm	14 in / 35.5 cm	350 lb / 158.8 kg	
7931KD-4, 7931KD-1	—	20 in / 50.8 cm	19 in / 48.2 cm	14 in - 20 in 35.5 cm - 50.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	14 in - 20 in 35.5 cm - 50.8 cm	—	350 lb / 158.8 kg	
7921R, 7921RB, 7921R-1, 7921RB-1	—	20 in / 50.8 cm	19 in / 48.2 cm	27 in - 33 in 68.5 cm - 83.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	14 in - 20 in 35.5 cm - 50.8 cm	14 in / 35.5 cm	350 lb / 158.8 kg	✓
7931R, 7931RB, 7931RB-1	—	20 in / 50.8 cm	19 in / 48.2 cm	14 in - 20 in 35.5 cm - 50.8 cm	20 in / 50.8 cm	12 in / 30.4 cm	14 in - 20 in 35.5 cm - 50.8 cm	—	350 lb / 158.8 kg	✓
7921RA	21 in / 53.3 cm	24 in / 60.9 cm	18 in / 45.7 cm	31 in / 78.7 cm	20 in / 50.8 cm	12 in / 30.4 cm	17 in - 20 in 43.1 cm - 50.8 cm	13 in / 33 cm	350 lb / 158.8 kg	✓
7931RA	21 in / 53.3 cm	24 in / 60.9 cm	18 in / 45.7 cm	25.2 in / 64.1 cm	20 in / 50.8 cm	12 in / 30.4 cm	17 in - 20 in 43.1 cm - 50.8 cm	—	350 lb / 158.8 kg	✓
7921R-BSKT-1	—	12 in / 30.4 cm	4 in / 10.1 cm	4 in / 10.1 cm	—	—	—	—	—	✓

Specifications are subject to change without notice

* EVENLY DISTRIBUTED

Bariatric

Lumex® provides a complete line of bathroom safety products specifically designed to meet the unique demands and weight limits of bariatric patients.

Imperial Three in One Drop Arm Steel Commode

- Extra-wide platform seat is comfortable for long periods of sitting
- Easy to release drop arm that locks automatically in raised position
- Commode opening accommodates the needs of larger users
- Item No. 6438A



Imperial Three in One Steel Commode

- Extra-wide platform seat is comfortable for long periods of sitting
- Commode opening accommodates the needs of larger users
- Item No. 7446A-2



Three in One Folding Steel Commode

- Folds for convenient storage and transport
- Sturdy plastic seat and lid snap on/off easily for cleaning
- Item No. 7109A-2



Three in One Aluminum Commode

- Aluminum frame is lightweight and easy to transport
- Contoured armrests enhance user comfort
- Sturdy plastic seat and lid snap on/off easily for cleaning
- Item No. 2190A



Transfer Bench

- Lightweight, durable, rust-resistant anodized aluminum frame
- Extra-deep seat provides added stability
- Drainage holes in seat minimize water build-up
- Item No. 7925A



Bath Seats

- Lightweight, durable, rust-resistant anodized aluminum frame
- Seat height adjustable in 1" increments to accommodate patient's size
- Available with or without backrest
- Item No. 7939A (With Backrest)
- Item No. 7932A (Without Backrest)



Item No.	Width Between Arms	Overall Width	Overall Depth	Backrest Height (Above Seat)	Seat Width	Seat Depth	Seat Height	Max. Weight Capacity *
6438A	25 in / 63.5 cm	31 in / 78.7 cm	22 in / 55.8 cm	-	23 in / 58.4 cm	18 in / 45.7 cm	19.5 in - 23.5 in 49.5 cm - 59.6 cm	600 lb / 272.1 kg
7446A-2	23 in / 58.4 cm	26 in / 66 cm	23.5 in / 59.6 cm	-	23 in / 58.4 cm	18 in / 45.7 cm	16.5 in - 23.5 in 41.9 cm - 59.6 cm	400 lb / 181.4 kg
7109A-2	24.5 in / 62.2 cm	30 in / 76.2 cm	20.5 in / 52 cm	-	14 in / 35.5 cm	16.5 in / 41.9 cm	17.5 in - 21.5 in 44.4 cm - 54.6 cm	650 lb / 294.8 kg
2190A	22.2 in / 56.5 cm	27 in / 68.5 cm	19 in / 48.2 cm	-	14.2 in / 36.1 cm	16.7 in / 42.5 cm	18 in - 23 in 45.7 cm - 58.4 cm	400 lb / 181.4 kg
7925A, 7925A-1	-	29.5 in / 74.9 cm	20.5 in / 52 cm	14.2 in / 36.1 cm	26.7 in / 67.9 cm	16.2 in / 41.2 cm	16.5 in - 20.5 in 41.9 cm - 52 cm	600 lb / 272.1 kg
7939A, 7939A-1	-	18 in / 45.7 cm	20.5 in / 52 cm	13 in / 33 cm	17 in / 43.1 cm	16.2 in / 41.2 cm	16 in - 20 in 40.6 cm - 50.8 cm	600 lb / 272.1 kg
7932A, 7932A-1	-	18 in / 45.7 cm	16.7 in / 42.5 cm	-	17.7 in / 45 cm	12 in / 30.4 cm	14 in - 17 in 35.5 cm - 43.1 cm	600 lb / 272.1 kg

Specifications are subject to change without notice

*EVENLY DISTRIBUTED

Hand-Held Shower Heads

The Lumex® Hand-Held Shower Heads provide the ultimate bathing experience. The simple controls make changing spray settings easy using just your fingertips even with wet soapy hands. Bathing pets, children, and persons with special needs is made easy with multiple spray options including a pause spray function. The Lumex Hand-Held Shower Head is easy to install on any standard shower arm within minutes. These hand-held shower heads come in attractive retail packaging.


Deluxe Hand Held Shower Head

- 6 spray options: saturating, saturating and bubbling, bubbling, massage, saturating and massage, and pause spray function
- 3-Way diverter valve
- Matching 84" reinforced hose provides extra convenience
- Item No. 12041 (Chrome finish)
- Item No. 12042 (Brushed nickel finish with stainless steel hose)



Everyday Hand Held Shower Head

- 5 spray options: saturating, bubbling, massage, half saturating, and pause spray function
- Includes mounting bracket and 2-way diverter valve
- Matching 80" reinforced hose provides extra convenience
- Item No. 12040



Universal Hand Held Shower Head Holder

- Fits almost any hand held shower head
- Angle-adjustable for easy reach
- Dual Mounting options: Drill-free suction cup or Screw-in
- Item No. 12040-HLDR-1 1 ea - Retail Packaging

Tub Rails

Lumex® Tub Rails provide a stable gripping surface for exit and entry into a bathtub. These white colored rails hold firmly to tub walls with clamps and will not mar the tub surface. They come in an attractive retail package.

Sure-Safe® Bathtub Rail

- Heavy-gauge steel is powder coated for rust resistance
- Provides two gripping surfaces to facilitate both entry and exit from tub
- Tool-free installation
- Fits tubs with 2.8" - 6.5" thick walls
- 250 lb maximum weight capacity, EVENLY DISTRIBUTED
- Item No. 2032A 1 ea - Retail Packaging



Heavy Duty Bathtub Safety Rail

- Textured, contoured handle provides a good no-slip surface
- Tall model has two gripping surfaces to facilitate both entry and exit from tub
- Tool-free installation
- Fits tubs with 2.5" - 6.3" thick walls
- 250 lb maximum weight capacity, EVENLY DISTRIBUTED
- Item No. 69623A (Standard Height - 12") 3/cs - Retail Pckg
- Item No. 69663A (Tall Height - 16.5") 3/cs - Retail Packaging



Transfer Benches

The Lumex® Transfer Bench is an excellent way to provide assistance to those who have difficulty stepping over a bathtub wall. These benches are lightweight, durable, and have rust-resistant aluminum frames. The reversible backrest accommodates any bathroom. Side support arm aids in lateral transfer.

Maxi-Drain Transfer Bench

- Bath seat has drain holes to minimize water build up
- Seat height is adjustable in 0.5" increments
- Item No. 7927A, 7927A-1



Knock-Down Transfer Bench

- Tool-free design for quick setup
- Designed for quick storage and easy transport
- Bath seat has drain holes to minimize water build up
- Seat height is adjustable in 0.5" increments
- Item No. 7927KD-2, 7927KD-1



Imperial Transfer Bench

- Larger and deeper seat than your standard transfer bench
- Durable blow-molded plastic bench and backrest
- Seat has curtain cut-out and accessory tray
- Item No. 7929, 7929-1



Padded Transfer Benches

- Comfortable, cushioned seat and backrest are watertight and easy to clean
- Seat height is adjustable in 1" increments
- Comes with tool-free assembly and back attachment
- Item No. 7955KD-1, 7955KD-2
- Item No. 7956KD-1 (with commode seat)



Item No.	Overall Width	Overall Depth	Seat Width	Seat Depth	Overall Height	Seat Height Range	Clearance Beneath Frame	Back Support (Above Seat)	Maximum Weight Capacity*
7927A, 7927A-1	29.2 in / 74.2 cm	22 in / 55.8 cm	26.5 in / 67.3 cm	16 in / 40.6 cm	31.7 in - 35.7 in 80.6 cm - 90.8 cm	18.2 in - 22.2 in 46.3 cm - 56.5 cm	15.2 in - 19.2 in 38.7 cm - 48.8 cm	13.5 in / 34.2 cm	350 lb / 158.8 kg
7927KD-2, 7927KD-1	29.5 in / 74.9 cm	17 in / 43.1 cm	25.5 in / 64.7 cm	16 in / 40.6 cm	33.5 in - 37.5 in 85 cm - 95.2 cm	18.2 in - 22.2 in 46.3 cm - 56.5 cm	16.2 in - 20.2 in 41.2 cm - 51.4 cm	14.5 in / 36.8 cm	350 lb / 158.8 kg
7929, 7929-1	32.5 in / 82.5 cm	23 in / 58.4 cm	28 in / 71.1 cm	16.7 in / 42.5 cm	30.5 in - 34.5 in 77.4 cm - 86.7 cm	17.5 in - 21.5 in 44.4 cm - 54.6 cm	15 in - 19 in 38.1 cm - 48.2 cm	13 in / 33 cm	400 lb / 181.4 kg
7955KD-1, 7955KD-2	29 in / 73.6 cm	22.7 in / 57.7 cm	24 in / 60.9 cm	15.7 in / 40 cm	32 in - 36 in 81.2 cm - 91.4 cm	17.7 in - 21.7 in 45 cm - 55.2 cm	14.5 in - 18.5 in 36.8 cm - 46.9 cm	13.7 in / 34.9 cm	300 lb / 136 kg
7956KD-1	29 in / 73.6 cm	22.7 in / 57.7 cm	24 in / 60.9 cm	15.7 in / 40 cm	32 in - 36 in 81.2 cm - 91.4 cm	17.7 in - 21.7 in 45 cm - 55.2 cm	14.5 in - 18.5 in 36.8 cm - 46.9 cm	13.7 in / 34.9 cm	300 lb / 136 kg

Specifications are subject to change without notice

*EVENLY DISTRIBUTED

Raised Toilet Seats & Risers

The Lumex® Raised Toilet Seats are designed to provide comfort and stability to those who have difficulty sitting down or standing up from a toilet. These seats, are designed to fit most toilets.

Everyday Raised Toilet Seat

- Lightweight and portable
- Adds 4.5" to toilet seat height
- Item No. 6909A, 6909A-1



Locking Raised Toilet Seat

- Front clamp locks securely to the toilet to prevent seat from shifting
- Wide, contoured seating surface is comfortable
- Available with or without removable padded armrests
- Adds 4.5" to toilet seat height
- Item No. 6486A, 6486A-1, 6486R (Without Armrests)
- Item No. 6487RA, 6487RA-1 (With Removable Armrests)



Multi-Position Open Padded Raised Toilet Seat

- Open seat promotes patient independence and facilitates perineal cleaning
- Front opening is 4.8" wide
- Seat can be rotated or rear mounted for maximum access to perineal area
- Thick, padded heat-sealed seat is ideal for patients with circulatory problems or who need additional comfort
- Comes with four heavy-duty VersaGuard™ coated locking brackets
- Item No. 6497A



Deluxe Toilet Seat Riser

- Adds 3.5" to toilet seat height
- Padded, removable armrests provide comfort
- Securely attaches to toilet bowl with screw and nuts
- Designed to function with your existing toilet seat to blend in with décor
- Designed for round or elongated toilets
- Item No. 6482R-2 (For Round Toilets)
- Item No. 6482ER-2 (For Elongated Toilets)



Toilet Safety Rails

- Helps support and enhance safety while toileting
- Armrests provide comfort and an excellent gripping surface
- Arms width adjusts for proper fit
- Item No. 6450R-2
- Item No. 6460A, 6460R, 6465A-1



Item No.	Width Between Arms	Overall Width	Overall Depth	Seat Width	Seat Depth	Seat Height Range	Weight Capacity*	Retail Pckg
6909A, 6909A-1	-	15 in / 38.1 cm			16.5 in / 41.9 cm	4.5 in / 11.4 cm	250 lb / 113.3 kg	
6486A, 6486A-1	-	16.5 in / 41.9 cm	17 in / 43.1 cm		16.5 in / 41.9 cm	4.5 in / 11.4 cm	300 lb / 136 kg	
6487RA, 6487RA-1	18.5 in / 46.9 cm	20.7 in / 52.7 cm	17 in / 43.1 cm		16.2 in / 41.2 cm	4.5 in / 11.4 cm	300 lb / 136 kg	✓
6497A	-	16 in / 40.6 cm	15 in / 38.1 cm	16 in / 40.6 cm	15 in / 38.1 cm	6 in - 7.5 in 15.2 cm - 19 cm	300 lb / 136 kg	
6060R, 6450R-2	18 in - 21 in 45.7 cm - 53.3 cm	20 in - 23 in 50.8 cm - 58.4 cm	16.2 in / 41.2 cm	-	-	-	250 250 lb / 113.3 kg	✓
6460A, 6465A-1	18 in - 21 in 45.7 cm - 53.3 cm	20 in - 23 in 50.8 cm - 58.4 cm	16.2 in / 41.2 cm	-	-	-	250 lb 250 lb / 113.3 kg	
6482R-2	19.5 in / 49.5 cm	22.5 in / 57.1 cm	17.2 in / 43.8 cm	13.7 in / 34.9 cm	17.2 in / 43.8 cm	3.5 in / 8.8 cm	300 lb / 136 kg	✓
6482ER-2	20.5 in / 52 cm	19.5 in / 49.5 cm	19 in / 48.2 cm	13.7 in / 34.9 cm	19 in / 48.2 cm	3.5 in / 8.8 cm	300 lb / 136 kg	✓

Specifications are subject to change without notice

* EVENLY DISTRIBUTED

Shower Transport Chairs

The Lumex® Shower Transport Chairs are designed for safe transport to and from a shower. These chairs roll easily over most toilets and include a comfortable nylon backrest that is removable for easy cleaning. The transport chairs are designed to be lightweight, durable, and have a rust-resistant frame.



gf
Proudly Made in the USA

PVC Shower Transport Commode Chair

- Constructed of the highest quality grade PVC
- Deluxe elongated open-front flip-up seat for easy cleaning
- Highest quality casters are made to roll easily over surfaces
- Quick drying back support is fluid-resistant, anti-bacterial, and anti-fungal
- Item No. 89180, 89220

Deluxe Shower Transport Chair

- Lightweight, durable, rust-resistant anodized aluminum frame
- Available with or without footrests
- Footplates tip up and out of away for easy transfer
- Item No. 7910A-1, 7915A-1




ITEM No.	Width between arms	Overall Width	Overall Depth	Overall Height	Seat Width	Seat Depth	Seat Height	Max. Weight Capacity*
89180	18 in / 45.7 cm	21.5 in / 54.6 cm	21.5 in / 54.6 cm	41.2 in / 104.7 cm	14.5 in / 36.8 cm	16.5 in / 41.9 cm	19.5 in / 49.5 cm	300 lb / 136 kg
89220	22 in / 55.8 cm	21.5 in / 54.6 cm	21.5 in / 54.6 cm	41.2 in / 104.7 cm	14.5 in / 36.8 cm	16.5 in / 41.9 cm	19.5 in / 49.5 cm	375 lb / 170 kg
7910A-1	17.5 in / 44.4 cm	21.5 in / 54.6 cm	24 in / 60.9 cm	39 in / 99 cm	14 in / 35.5 cm	16.5 in / 41.9 cm	22 in / 55.8 cm	250 lb / 113.3 kg
7915A-1	17.5 in / 44.4 cm	21.5 in / 54.6 cm	33 in / 83.8 cm	39 in / 99 cm	14 in / 35.5 cm	16.5 in / 41.9 cm	22 in / 55.8 cm	250 lb / 113.3 kg

Specifications are subject to change without notice

* EVENLY DISTRIBUTED

Grab Bars

Lumex® Grab Bars are constructed of heavy duty steel with a rust-resistant finish. They are available in a chrome plated, or VersaGuard™ (white enamel) finish. A knurled surface provides a good gripping surface and they can be mounted vertically, horizontally, or diagonally. These grab bars carry a 250 lb evenly distributed maximum weight capacity. They come in an attractive retail package.



VersaGuard® Coated

Item No.	Length	Qty
40033A	12 in / 30.4 cm	3/cs
40433A	16 in / 40.6 cm	3/cs
40133A	18 in / 45.7 cm	3/cs
40233A	24 in / 60.9 cm	3/cs
40333A	32 in / 81.2 cm	3/cs



Chrome Plated

Item No.	Length	Qty
2012A	12 in / 30.4 cm	1ea
2016A	16 in / 40.6 cm	1ea
2018A	18 in / 45.7 cm	1ea
2019A	24 in / 60.9 cm	1ea
2021A	32 in / 81.2 cm	1ea



Follow us on

1.770.368.4700

Information contained herein is subject to change. The most current and complete product information can be found on our website.
www.grahamfield.com

© 2013, GF Health Products, Inc. All Rights Reserved. Lumex is a trademark of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.