



Lifts and Slings
Wheelchairs
Bath Safety





A Graham-Field® Family Brand

A safe lifting program is essential in today's healthcare environment. Our population is getting older, heavier, and has more chronic health conditions than ever before. This puts an increasing demand on our caregivers. Using lifts is safer for patients and healthcare personnel.

Lumex lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use. All Lumex floor lifts have a manual lowering mechanism.

PRO Battery-Powered Floor Lift



Premium electronics with LCD Display



Control box AND intuitive hand pendant operate high/low and base open/close



Battery charge indicator on hand pendant and LCD Display



Includes charging caddy, charger, and spare battery



Alerts clearly shown on LCD Display



Swan neck enables lifting resident from the floor



Optional digital scale (DSC270)

Battery-Powered Floor Lifts

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Large, easy to grip handles for easy maneuverability
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)



Ergonomic foot pedal leg opening frees up hands



Charger break-away cord



Storage bag included



Optional digital scale (DSC250)

Battery-Powered Sit-to-Stand Lifts

Sit-to-stand lifts are ideal for patient transfers when the patient is actively and mentally able to participate and is weight bearing in both legs.

They are designed for patients who have difficulty standing up on their own from a seated position and the reverse.



The Lumex Stand Assist is a transport assistance unit that keeps the resident actively engaged in the transport process. Transport functions of all kinds are quick, convenient, and require minimal caregiver assistance. A padded split seat swings out to allow loading or unloading, then swings back to form a

comfortable and stable transport seat. The Stand Assist is a convenient toileting alternative to a wheelchair. *The Lumex Stand Assist is designed for residents who have the strength and stability to lift and support themselves. We recommend a battery-powered sit-to-stand lift for residents with limited strength and stability.*



The user grasps the middle bar and uses their own strength to pull themselves up into position.





Full Body Slings

Full Body slings are the easiest sling to use while transferring a resident from a supine position (laying on ones' back). It can be used for transfers from a bed to a chair, wheelchair, or toilet, or from the floor to a bed. This sling style offers total head, neck, and thigh support. Full Body slings should be used for residents who are totally or partially dependent, non-weight bearing, bariatric, or for those who have limited head control. These

are available in a solid fabric, mesh, or mesh with a commode opening. Mesh slings should be used for bathing. These slings are often left under the resident between lifts. Solid fabric slings should not be left under the resident for extended periods of time, as the fabric can retain body heat. LUMEX Full Body slings have four connection points and can be used with 2-, 4-, or 6-point spreader bars.



Solid Fabric



Mesh



Mesh with Commode Opening

Universal Slings

Universal slings (also called divided leg slings) are the easiest sling to use while transferring a resident from a seated position. These can also be used for lifting from a supine position, making them suitable for most lifting situations. They are an easy-to-fit general purpose sling designed for the majority of users. They can be used for transfers from a bed to a chair, wheelchair, or toilet, chair to chair, or from the floor to a bed or chair. This sling style is available with or without head support. They are available in solid fabric and mesh styles.

Universal slings should be used for residents that are totally or partially dependent, non-weight bearing, or bariatric. Those who have limited head control should use a universal sling with head support. The large opening accommodates toileting. Mesh slings should be used for bathing. LUMEX Universal slings have four or six connection points and can be used with 2-, 4-, or 6-point spreader bars.



Solid Fabric



Solid Fabric with Head Support



Mesh



Mesh with Head Support

Toileting Slings

Toileting slings are designed to facilitate the toileting procedure by allowing the lowering of clothing. The sling is well padded under the arms and in the lumbar area for maximum comfort. A hook and loop belt and separate buckle strap hold the sling securely in place. Toileting slings must only be used when transferring a resident from a seated position. They can be used for transfers

from sitting on a bed to a chair, wheelchair, or toilet, or chair to chair.

LUMEX Toileting slings have four connection points and can be used with a Sit-to-Stand lift or 2-, 4-, or 6-point spreader bars.



Deluxe Sit-to-Stand Padded Slings

Sit-to-Stand slings are designed to facilitate transfers using sit-to-stand lifts. It is only suitable for residents who have a greater degree of weight bearing ability. In addition to making seating transfers, this sling will facilitate the toileting procedure by allowing the lowering of clothing. The sling is well padded under the arms and in the lumbar area for maximum comfort. A textured belt and separate buckle strap hold the sling securely in place. Sit-to-Stand slings must only be used when

transferring a resident from a seated position. They can be used for transfers from sitting on a bed to a chair, wheelchair, or toilet, or chair to chair. LUMEX Sit-to-Stand slings have two connection points and can be used only with a Sit-to-Stand lift.

Sit-to-Stand slings can also be used in conjunction with an optional Buttock Strap and/or a Knee/Calf Strap for extra security and support.



Wheelchairs



Everest & Jennings®

A Graham-Field® Family Brand

Has a long tradition of quality wheelchairs, a brand respected for performance and durability.

We offer a full range of wheelchairs to meet your needs.



Traveler SE
Standard Hemi



Traveler LX
Premium Chrome-plated



Traveler L3
Lightweight



Traveler L4
High-strength
Lightweight



Traveler HD
Heavy Duty



**Advantage
Recliner**



Paramount XD
Extra Heavy Duty

Bath Safety

LUMEX®

A Graham-Field® Family Brand

We offer a full range of our LUMEX bath safety products for your facility.



Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame



Bath Seats with and without Arms



Bariatric Bath Seats



Transfer Benches



PVC Shower Chairs/
Shower Gurney

Lifts and Slings Specifications

Model	Floor Clearance	Spreader Bar Minimum Height from Sling Attachment	Spreader Bar Maximum Height from Sling Attachment	Overall Length	Minimum Base Width (outside)	Maximum Base Width (outside)	Safe Working Load (EVENLY DISTRIBUTED)	Product Weight
LF500	4.3"	19.3"	74.8"	49.2"	25.6"	39.4"	500 lb	116.8 lb
LF1050	4.9"	27.5"	68.1"	40.2"	26.0"	37.8"	400 lb	114.5 lb
LF1090	4.9"	27.2"	70.9"	50.0"	26.0"	43.0"	600 lb	137.0 lb
LF2020	4.9"	36.4"	70.8"	40.4"	25.5"	37.8"	400 lb	118.0 lb
LF2090	4.9"	36.4"	70.8"	40.4"	25.5"	37.8"	600 lb	134.0 lb

Model	Safe working load	User Height	Maximum Width	Maximum Height	Overall Length	Unit Weight	Front Base Height
LF1600	400 lb	Between approximately 5' 1" and 6' 6"	25"	42.5"	33.5"	61.0 lb	4.7"

Full-Body Slings	Style	Item	Size	Accent Color	Safe Working Load (lbs.)	Best Fits (lbs.)
	Solid Fabric	F112	M	Red	450	100-195
F113		L	Green	450	150-300	
F117		XL	Blue	450	200-400	
Mesh	FM110	M	Red	450	100-195	
	FM111	L	Green	450	150-300	
	FM140	XXL	Yellow	600	290-600	
Mesh with Commode Opening	FMC114	M	Red	450	100-195	
	FMC115	L	Green	450	150-300	
	FMC116	XL	Blue	450	200-400	
	FMC141	XXL	Yellow	600	290-600	

Universal Slings	Style	Item	Size	Accent Color	Safe Working Load (lbs.)	Best Fits (lbs.)
	Solid Fabric	UP892	M	Yellow	400	140-200
UP893		L	Green	400	200-400	
UP894		XXL	White	600	400-600	
Solid Fabric with Head Support	UPH862	M	Yellow	400	140-200	
	UPH863	L	Green	400	200-400	
	UPH864	XXL	White	600	400-600	
Mesh	MS882	M	Yellow	400	140-200	
	MS883	L	Green	400	200-400	
Mesh with Head Support	UMH832	M	Yellow	400	140-200	
	UMH833	L	Green	400	200-400	

Toileting Slings	Style	Item	Size	Accent Color	Safe Working Load (lbs.)	Torso (inches)
	Deluxe Padded	TP332	M	Yellow	400	30-48
TP333		L	Green	400	38-60	
TP334		XXL	White	600	58-73	

Sit-to-Stand Slings	Style	Item	Size	Accent Color	Safe Working Load (lbs.)	Torso (inches)
	Deluxe Padded	DLSA1	S	Red	400	26-42
DLSA2		M	Yellow	400	30-46	
DLSA3		L	Green	400	34-53	
DLSA4		XL	White	600	46-64	
Buttock Strap	DLSA9	Std	N/A	400	N/A	
	DLSA9B	Bari	N/A	600	N/A	
Knee/Calf Strap	SAL1230-KB	Std	N/A	N/A	N/A	



1.770.368.4700



Information contained herein is subject to change. The most current and complete product information can be found on our website.

www.grahamfield.com

© 2017, GF Health Products, Inc. All Rights Reserved.

Lumex, E&J, Everest & Jennings, Graham-Field, and Traveler are trademarks of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.