

DISTRIBUTED BY:



# TrueClr M<sup>+</sup>

## Active External Urine Collection Apparatus

---

TrueClr - Preserve Dignity & Save Time

# AN ALTERNATIVE TO USING A TRADITIONAL INTERNAL CATHETER

TrueClr is the world's first completely external system designed to comfortably evacuate the bladder while eliminating CAUTIs.

- Non-invasive & non-adhesive Eliminates skin irritation & inflammation
- Latex-free & hypoallergenic
- Easy-to-use, portable & reusable
- Made in the USA

## WHAT YOU NEED

1. A TrueClr Urine Collection Apparatus (UCA)
2. An aspirator pump. Aspirator must meet specifications:
  - Vacuum Range 50-550 mm Hg +/- 10%\*
  - Air Flow @ Pump Inlet: 27 LPM (free flow)
  - Reusable Collection Container Capacity 1200 ml (cc)


**Figure 1**

Roll back the end of UCA



**Figure 2**

Apply the UCA and unroll onto penis



**Figure 3**

Urine flows into canister



**Figure 4**

Break suction and unroll UCA to remove



## TrueClr M+

A revolutionary external catheter system designed to actively empty the bladder while eliminating the risk of Catheter Associated Urinary Tract Infections (CAUTIs)

**>13,000**

deaths per year from internal catheters

**80%**

of all hospital Urinary Tract Infections (UTIs) result from internal catheters

**\$450**

spent in healthcare due to CAUTIs in the USA



### Patents

62/297,096 008275770-0002 15/416,272 2020-025534 11,246,573 1706212 35/201,036 2021016163  
15/412,049 008275770-0003 US 2017/0238911 1713030 17/586,055 29/736,239 US 11,311,405 202016166  
US 2017/0241978 008275770-0004 US 10,682,124 2020-25535 US 2022/0202398 29/825,375 17/657,658 199397  
US 10,690,655 2020-025533 16/901,614 1713031 JP PAT # 6,895,182 202016158 16/681,601 PCT/AU2021/020930 212723  
JP PAT # 6,839,452 1713029 US 2020/0345332 2020-25536 202016162 WO 2021/236202 008275770-0001