

BOSMA
Enterprises®

OPPORTUNITY PRO™

NITRILE EXAM GLOVES

with Compression Fit Technology

Opportunity Pro nitrile gloves from Bosma Enterprises are **manufactured in the USA** with TAA compliant foreign NBR, and designed for extended wear and use throughout a facility, from exam rooms to environmental management. Opportunity Pro gloves are comfortable while providing protection from **fentanyl** and chemotherapy drugs. These gloves include **Compression Fit Technology**, a patented manufacturing technology that produces compression in fingertips for improved comfort, enhanced dexterity, and greater tactile sensitivity, without compromising durability or protection.

Packaging: 100/box; 10 boxes/case (XS, S, M, L and XL)
90/box; 10 boxes/case (2XL)

Color: Blue

Finish: Textured fingertip

Type: Ambidextrous/Non-Sterile/Disposable/Single Use

Length: 9.5 inches

Cuff: Beaded

Palm thickness: 0.08 mm

Finger thickness: 0.135 mm

Tensile Strength: >35 MPa

Elasticity/elongation: >550%

Quality Standards:

- Tested for use with chemotherapy drugs in accordance with ASTM D6978
- Barrier tested against fentanyl exposure per ASTM D6978-05
- 1.5 AQL
- ISO 9001 and 13485
- 510(K) / Approved for medical use
- Component materials comply with FDA regulations for food contact

Size

X-Small	72-50800
Small	72-50801
Medium	72-50802
Large	72-50803
X-Large	72-50804
2X-Large	72-50805



Manufactured by American Nitrile
Grove City, OH



For more information
go to bosma.org or
call **800.362.5463**

BOSMA
Enterprises®

6270 Corporate Drive
Indianapolis, IN 46278

Bosma Enterprises®, creating opportunities for individuals who are blind or visually impaired.
Your purchase creates employment for more than 100 people who are blind.



ISO 9001:2015
ISO 13485:2016
CERTIFIED

