

# BOSMA Enterprises®

## OPPORTUNITY FLEX™

NITRILE EXAM GLOVES

with Compression Fit Technology

Opportunity Flex nitrile gloves from Bosma Enterprises are **manufactured in the USA** with TAA compliant foreign NBR and designed for medical professionals that wear gloves all day with frequent changes. Opportunity Flex gloves offer the perfect blend of softness, strength, high tactile sensitivity, and easy donning, while providing protection from **fentanyl** and chemotherapy drugs. These gloves include **Compression Fit Technology**, a patented manufacturing technology that produces compression in fingertips for improved comfort, enhanced dexterity, and greater tactile sensitivity, without compromising durability or protection.

**Packaging:** 200/box; 10 boxes/case (XS, S, M and L)  
180/box; 10 boxes/case (XL and 2XL)

**Color:** Periwinkle Blue

**Finish:** Textured fingertip

**Type:** Ambidextrous/Non-Sterile/Disposable/Single Use

**Length:** 9 inches

**Cuff:** Beaded

**Palm thickness:** 0.07 mm

**Finger thickness:** 0.11 mm

**Tensile Strength:** >30 MPa before aging; >27 MPa after aging

**Elasticity/elongation:** >500% before aging; >400% after aging

### Quality Standards:

- Tested for use with chemotherapy drugs in accordance with ASTM D6978
- Barrier tested against fentanyl exposure per ASTM D6978-05
- 1.5 AQL
- ISO 9001 and 13485
- 510(K) / Approved for medical use
- Component materials comply with FDA regulations for food contact

### Size

X-Small	72-50900
Small	72-50901
Medium	72-50902
Large	72-50903
X-Large	72-50904
2X-Large	72-50905



Manufactured by American Nitrile  
Grove City, OH



For more information  
go to [bosma.org](http://bosma.org) or  
call **800.362.5463**

**BOSMA Enterprises®**

6270 Corporate Drive  
Indianapolis, IN 46278

Bosma Enterprises®, creating opportunities for individuals who are blind or visually impaired.  
Your purchase creates employment for more than 100 people who are blind.



ISO 9001:2015  
ISO 13485:2016  
CERTIFIED

